



Strip	
1	anti-MuLV (1:2000)
2	anti-MuLV (1:2000) + p15 (100 ug/ml)
3	anti-MuLV (1:2000) + p15E (100 ug/ml)
4	anti-MuLV (1:2000) + p30 (200 ug/ml)
5	anti-MuLV (1:2000) + p10 (100 ug/ml)
6	anti-MuLV (1:2000) + p12 (44 ug/ml)
7	anti-Env (1:1000)
8	anti-Env (1:1000) + gp70 (30 ug/ml)